

UPCOMING EVENTS

**Beyond Traditional Lipids:
Case-Driven Cardiometabolic
Risk Assessment with
Advanced Biomarkers**

Presented by Melody Hartzler, PharmD
March 11, 2026 at 12 PM Pacific

**Translating Methylation
Science into Actionable
Insights: Lessons from the
Doctor's Data Plasma
Methylation Profile**

Presented by Julia Malkowski, ND, DC
April 1, 2026 at 12 PM Pacific

Visit doctorsdata.com/Register-for-Webinars




WILL BEGIN SHORTLY




Hormone Metabolism Simplified: A Practical HuMap™ Guide for Busy Integrative Clinicians

Presented Ruth Hobson, ND



The Hormone and Urinary Metabolites Assessment Profile (HuMap™)

- Provides a comprehensive overview of steroid hormones, their metabolites, and the efficiency of the enzymes that metabolize these hormones.
- This non-invasive test requires only 4 or 5 separate liquid urine collections.
- As the liver and digestive organs play a role in hormone breakdown and elimination, urinary testing of hormones and metabolites may elucidate a need to address other areas of health.
- Additionally, testing urinary hormone metabolites can contribute to further understanding of endogenous hormone secretion, supplemental hormone utilization, enzyme activity, oxidative stress, and insight into how the body metabolizes hormones.




How are urinary hormones measured?

LC/MS/MS

- LC separates hormones in a liquid urine sample, which are then injected at various times for MS analysis.
- MS technology monitors the injected sample for specific hormones based on their molecular weights and expected injection times (commonly called retention times).
- Tandem MS verifies each hormone identity based on fragmentation and determines its amount.
- The combination of LC and tandem MS allows for extremely sensitive and specific hormone measurements, even in samples containing similar substances that would interfere with other methods.

What is measured in urinary hormone testing?



Unconjugated hormones:


Sex hormones like progesterone, testosterone, estradiol, etc.

Conjugated metabolites:

Consist of hormones that the body has modified by adding a water-soluble side arm - usually glucuronide or sulphate molecules - to facilitate their elimination in urine or stool.

What are the benefits of looking at hormone urine metabolites over salivary hormones?

Which has been said to be more accurate?



When to consider using HuMap™

General

- 🔗 Metabolic syndrome
- 🔗 Thyroid pathologies
- 🔗 Inflammation
- 🔗 Oxidative Stress
- 🔗 Fatigue/insomnia
- 🔗 Low libido
- 🔗 Mood or cognitive concerns
- 🔗 Family/personal history of hormone-driven cancers
- 🔗 HRT/BHRT utilization
- 🔗 Weight gain
- 🔗 Alopecia



Female

- 🔗 Menopausal symptoms
- 🔗 Breast Health
- 🔗 Endometriosis
- 🔗 PCOS
- 🔗 PMS/PMDD

Male

- 🔗 Erectile dysfunction
- 🔗 Prostate health
- 🔗 Lower urinary tract symptoms
- 🔗 Loss of muscle mass
- 🔗 Breast health

Common actionable findings on HuMap reports

Androgens

- 🔍 Excess aromatase or 5 alpha reductase activity depleting testosterone in males
- 🔍 Excess 5 alpha reduced androgens in females causing symptoms

Progesterones

- 🔍 Lack of abundance of progesterones in premenopausal females suggesting anovulation (or in males, impacting prostate health)

Oxidative Stress

- 🔍 Elevated 8-OHdG

Estrogens

- 🔍 Imbalance between 2-, 4-, and 16-hydroxy estrone (Phase I)
- 🔍 Impaired methylation
- 🔍 Insufficient phase III detoxification (constipation/dysbiosis) leading to hormone recirculation and excess levels

Corticoids

- 🔍 Elevations at waking suggestive of excess cortisol excretion while sleeping
- 🔍 Imbalance between cortisol and cortisone
- 🔍 Increased or decreased excretion of corticoids

5 COLLECTING

Day 1 Collection		Midsleep
	#1 Dinnertime: collect within 1-2 hours prior to eating	<i>Midsleep: collect if you wake to urinate while sleeping</i>
	#2 Bedtime: collect within 1 hour prior to sleep	

- Collect a midstream specimen of urine provided. If needed, you may reuse the
- Pipette or pour urine into the properly marked on the tube label. Any remaining discarded in the toilet. DO NOT urinate an acid preservative* that may irritate the
- Close tube and rock gently.
- Record your name and the date and time on the requisition form.
- Place the urine tube back into the original absorbent pad. Store the urine in the freezer.

#1 Dinnertime (orange): Collect the first urine sample within 10 minutes.

#2 Bedtime (blue): Collect the second urine sample within 10 minutes.

Midsleep (grey): Collect only if you wake to urinate during your sleep.

#3 Waking (pink): Upon waking and rising from bed, collect the waking urine sample within 10 minutes.

#4 Post-waking (green): Wait at least 2 hours and collect another urine sample within 2-3 hours after waking. If you are unable to urinate, drink a little fluid.

Complete all the request information on the patient survey. All urine samples should be stored in the freezer with the frozen ice pack for 24 hours before shipping.

*Acid preservative is approximately 25mg Oxalic acid (found in many household cleaning products.) In case of contact with skin, wash immediately with water.

Day 1 Collection		Midsleep	Day 2 Collection	
	#1 Dinnertime: collect within 1-2 hours prior to eating	<i>Midsleep: collect only if you wake to urinate while sleeping</i>		#3 Waking: collect within 10 mins of waking
	#2 Bedtime: collect within 1 hour prior to sleep			#4 Post-waking: collect after 2-4 hours

- Collect a midstream specimen of urine using the collection cup provided.
- Pipette or pour urine into the properly labeled tube to the "fill line" marked on the tube label. Any remaining urine in the collection cup can be discarded in the toilet. Rinse the collection cup with clean water and allow to air dry before the next urine collection. DO NOT urinate directly into the tube, which contains an acid preservative* that may irritate the skin.
- Close tube and rock gently. Record your name and the date and time of collection on the tube, and on the requisition form.
- Place the urine tube back into the original plastic bag containing the absorbent pad. Store the urine in the freezer with the frozen ice pack. Discard the urine collection cup.

#1 Dinnertime (orange): Collect the first urine sample prior to eating dinner. Try to wait 3-4 hours from your last urination before this first collection.

#2 Bedtime (blue): Collect the second urine sample prior to going to bed – this should be the last urination for Day 1.

Midsleep (grey): Collect only if you wake to urinate during your sleep.

#3 Waking (pink): Upon waking and rising from bed, collect the waking urine sample within 10 minutes.

#4 Post-waking (green): Wait at least 2 hours and collect another urine sample within 2-4 hours after waking. If you are unable to urinate, drink a little fluid.

Hormone and Urinary Metabolites Assessment Profile

URINE COLLECTION



INSTRUCTIONS

Read all instructions before collecting specimens.



Scan code for expanded instructions and FAQs



How does a 4-point liquid urine collection method compare to 24-hour collection?



4 (or 5) - point urine collection throughout a 24-hour period offers a few advantages over 24-hour collection.



Collection of 4 urine samples is easier than collecting all urine in a 24-hour period as up to 40% of patients who attempt 24-hour urine collection do not collect properly.




4-point collection allows measurement of the diurnal cortisol rhythm whereas 24-hour collection can only provide a measure of overall cortisol production.



Hormone values measured from a pool of the 4 collections correspond well with 24-hour hormone values.

What are the advantages of a liquid sample?

- 🌱 The main advantage of liquid urine collection is enhanced sensitivity, especially for low concentration metabolites. Dried urine must be reconstituted from the filter paper once the sample arrives. This reconstitution can lead to loss of polar steroid metabolites or creatinine for some patient samples.
- 🌱 With liquid urine, samples can be shipped after being frozen for 4-6 hours, can be processed faster, and concentrated further to enhance the detection of low-level analytes. Steroid hormones and metabolites are also quite stable in liquid urine if the correct preservative is used.
- 🌱 Ability to add on other urinary tests (NTs or Amino Acids)



BEST PRACTICES FOR SPECIMEN COLLECTION

HORMONE AND URINARY METABOLITES ASSESSMENT PROFILE

The following collection recommendations are specific to Doctor's Data. Adhering to these recommendations ensures test results correlate with established reference ranges. *Never discontinue prescription medications without first consulting your provider.*

version 4.22

Hormone Supplementation	Timing of Last Supplementation
Topical, IM/SQ Injections, subcutaneous pellets, transdermal patches	Continue hormones based on your regular schedule
Vaginal	Discontinue 72 hours prior and during sample collection as this may directly contaminate the urine
Sublingual (dissolved under the tongue), Oral	Discontinue 72 hours prior and during sample collection
Cortisol/Hydrocortisone	Discontinue 4-5 days prior to sample collection to test endogenous secretion
Additional Considerations	
Glucocorticoid supplementation: Consult with your health care provider if you are taking glucocorticoids. Certain medications such as asthma inhalers and anti-itch creams contain cortisol and can impact cortisol/cortisone levels.	
24 hours prior and during collection: avoid alcohol, caffeine, tobacco or nicotine-containing products and strenuous exercise.	
On day 1 of collection, do not drink more than 4 liters of fluids.	
If you are also testing Neurotransmitters from the same urine samples	
24 hours prior and during collection: avoid avocados, eggplant, tomatoes, bananas, melons, pineapple, grapefruit, plums, fruit juice, nuts, nut butters, wine, cheese, rice, and chocolate.	
Day of both collections: it is recommended to avoid all supplements and medications until after all samples have been collected (including those that regulate allergy, mood, sleep, pain and inflammation)	

Collection Schedule

Menopausal Status	When to Collect Samples
Premenopausal, regular cycles (28 days)	Days 19-23 (Mid luteal phase)
Premenopausal, regular cycles longer than 28 days	Count back 7-9 days from usual end of cycle, at a minimum of day 19 (Luteal phase is almost always 14-16 days long)
Perimenopausal, irregular cycles with ovulation pains	7 days after ovulation
Perimenopausal, irregular cycles, no ovulation pains	Test after day 14, and before day 1 of next cycle
Premenopausal, irregular cycle shorter than 14 days	Days 7-9
When there is no point of reference, e.g. 60 days	Collect samples; freeze; if no period within 2 days then mail the sample
Men and postmenopausal women	Anytime

Schedule	When to Collect Sample
Night shift workers	Shift workers should adapt the collection times in accordance with their own sleep / wake schedule.
Frequent travelers (especially across multiple time zones)	If possible, collect samples after 2 weeks at home on regular schedule.

Health Disclaimer: All information given about health conditions, treatments, products and dosages are not intended to be a substitute for professional medical advice, diagnosis or treatment.





Hormone & Urinary Metabolites Assessment Profile



Order: 999999-9999



Test: X999999-9999-1

Client #: 999999

Doctors Data Inc

123 Main St.

St. Charles, IL 60174 USA

Patient: Sample Patient

Id: 999999

Age: 61 DOB: 01/01/1960

Sex: Female

Body Mass Index (BMI): 25

Menopausal Status: Post-menopausal

Sample Collection Date/Time

Dinnertime 12/30/2022 19:20

Bedtime 12/30/2022 22:30


Waking 12/31/2022 07:00

2 Hr. Post Waking 12/31/2022 09:30

Collection Period Multipoint daily

Date Received 01/01/2023

Date Reported 01/02/2023



ESTROGENS

The bar graph represents the relationship of the catechol estrogens (2-OH-E1, 4-OH-E1, 16-OH-E1) to each other. The expected percentage for each is represented by the shaded area.


The pathway illustrates phase 1 and phase 2 metabolism of both E1 and E2. Phase 1 metabolites, also known as catechol estrogens, are active and can induce estrogenic actions. Phase 2 metabolism gives insight into a patient's ability to methylate, or potentially inactivate harmful metabolites.



2-OH: generally considered safest



4-OH: potential for DNA damage

16-OH: considered highly estrogenic



CORTICOIDS

11 β HSD2 is responsible for the conversion of cortisol to cortisone. Inhibition of this enzyme may lead to the amount of cortisol being greater than cortisone, while increased enzyme activity can lead to higher levels of cortisone in comparison to cortisol.



KEY RELATIONSHIPS

The graphs to the right represent metabolism preference by key enzymes, indicated by the arrow.

Metabolites in the 5-alpha pathway are more androgenic than their 5-beta counterparts and can be responsible for androgenic symptoms even when hormone levels appear normal.

Aromatase is an enzyme found in the greatest amounts in peripheral fat tissue which can increase estrogens in both males and females.

4-OH-E1 is considered unfavorable due to its carcinogenic potential within breast and prostatic tissue as a reactive metabolite. When methylated by COMT, this reactive metabolite becomes stable and can be removed from the body.

5-A REDUCTASE ACTIVITY




AROMATASE ACTIVITY



COMT/METHYLATION ACTIVITY



HuMap Colors




- WRI – Within Reference Interval - represented by bracket and stated ranges on report
- Dark Blue = Below RI
- Light Blue = WRI low
- Green = Optimal
- Yellow = WRI high
- Red = Above RI

Neighborhoods

🌱 The HuMap™ test can be divided into 4 neighborhoods

- 🌱 Progesterones
- 🌱 Corticoids
- 🌱 Androgens
- 🌱 Estrogens





Progesterone Neighborhood



- 🔗 Progesterone
- 🔗 5-alpha pregnanediol
- 🔗 5-beta pregnanediol
 - 🔗 Dominant metabolite
 - 🔗 Can correlate with ovulation
- 🔗 Allopregnanolone
- 🔗 5-pregnenetriol
 - 🔗 Elevations seen in PCOS, CAH, Cushing's
- 🔗 17-hydroxyprogesterone
 - 🔗 Elevations seen in PCOS, idiopathic hirsutism, CAH, and NCCAH, hyperinsulinemia, stress, and alcohol use
- 🔗 21-hydroxyprogesterone
 - 🔗 Elevations seen in Cushing's, DMII, CAH

Allopregnanolone



- 🌀 Crosses the BBB and binds the GABA receptor
- 🌀 May help evaluate oral progesterone metabolism



Is there any utility to testing a patient on oral contraceptives?



23yo F, OCP
(E2 + progestin)




29yo F, progestin IUD



Corticoids

Corticoids Neighborhood

- Graphical representation of each timed collection
- Cortisol vs. Cortisone
- Metabolized cortisol
- Cortisol vs. cortisone metabolites




Diurnal Pattern



- 🕒 Time stamps correspond to time 3-4 hrs before collection
- 🕒 Gives information on overnight cortisol secretion.
 - 🕒 Can identify one etiology of insomnia/nighttime wakefulness
 - 🕒 Often associated with hypoglycemia overnight

Urinary Corticoids = HPA axis + Bodies Response to Stress

- 🌱 Cortisol (active) vs. Cortisone (storage)
 - 🌱 11BHS2 enzyme: most active in the kidney
 - 🌱 11BHS1-pulls cortisol out of storage
 - 🌱 Affected by obesity, metabolic syndrome, inflammation, thyroid function.
 - 🌱 11BHS2- protects the mineralocorticoid receptor from cortisol effects
- 🌱 Because we are looking at the influence of the kidney, we aren't necessarily seeing the HPA axis activity directly, rather the bodies response to cortisol levels



Corticoid Ratios and Calculations


- 🌱 Metabolized cortisol: (THE + 5 α THF + 5 β THF)
 - 🌱 Measure of what the body has utilized
 - 🌱 High: increased cortisol clearance due to over production of cortisol, also influenced by; obesity, IR, inflammation, hyperthyroidism
 - 🌱 Low: decreased cortisol clearance which could also be due to obesity, inflammation, IR, hypothyroidism
- 🌱 Cortisol metabolites/ Cortisone metabolites: (5 α THF + 5 β THF) vs. THE
 - 🌱 Metabolite ratios can give clinicians a better idea of overall preference within the body

Patterns seen on corticoid reports

Underweight patient may have decreased cortisol metabolites




Hypothyroidism can slow cortisol clearance and decrease cortisol metabolites



Obesity can increase cortisol clearance and 11βHSD1 activity causing a cortisol preference




Hyperthyroidism can upregulate liver clearance of cortisol



What downregulates 11B-HSD1 activity?

- Reduce inflammation, insulin resistance/insulin, central adiposity
- Physical activity
- Rooibos
- Holy Basil
- Curcumin
- Bitter melon
- EGCG
- Progesterone
- Coffee





Androgens

Can You Explain Aromatization?



Aromatase (CYP19A1)


- Converts androstenedione to estrone and testosterone to estradiol
- Active in adipose tissue
- Many of these CYP450 genes are subject to genetic polymorphisms
- Well-known for its role in the carcinogenic bioactivation of polycyclic aromatic hydrocarbons (PAHs), heterocyclic aromatic amines/amides, polychlorinated biphenyls (PCBs), and other environmental toxins




Activation of Aromatase (CYP19A1)

- Stress/cortisol: PMID: 21878510 PMID: 23835908 (animal studies, acts in a protective manner for brain health)
- Xeno-estrogens/endocrine disruptors/toxins (pesticides, herbicides, benzene, plastic by-products, some pharmaceuticals and cosmetics, petroleum, UV filters) PMID: 2223368
- Estrogen PMID: 2223368
- Poor dietary choices/high glycemic foods PMID: 2223368
- Major illness PMID 16670151
- Brain injury PMID: 16498364 (animal study)
- Alcohol/Red wine PMID: 19268535
- Forskolin (found in coleus plant) PMID: 14709151
- High insulin: PMID: 3322018
- Excess adipose tissue/Obesity/Leptin resistance: PMID: 10349800 PMID: 3226011
- Free-fatty acids (increased in metabolic syndrome): PMID: 2223368
- Inflammatory cytokines: IL-6, TNFalpha, prostaglandin PGE(2) PMID: 10405348

Aromatase activity and its effect on estrogen levels in males



63yo M with Metabolic Syndrome on transdermal testosterone causing gynecomastia!



Aromatase and Breast Cancer


- Intracellular androgen receptor and the membrane androgen receptor both down-regulate the strongly antiapoptotic protein Bcl-2. *protect
- Estrogen receptor- α and membrane estrogen receptor both up-regulate Bcl-2. *proliferate
- **All of this is consistent with the increase in breast cancer observed in postmenopausal women who had higher serum testosterone levels being due to the resulting increase in the local level of estradiol in the breast tissue. *proliferate**

Aromatase (CYP19A1)

CYP19A1


(-) Chrysin, Zn, nettles, damiana, grape seed, ECGC, resveratrol, licorice, flavonoids, anti-inflammatories

(+) Coleus forskohlii, inflammation, stress, Zn deficiency, excess adipose, high insulin, alcohol



5 α -reductase

- Enzyme that converts testosterone to its more potent form, dihydrotestosterone (5 α -DHT)
 - DHT has 2-3 times the affinity for the androgen receptor than testosterone
- May lead to increase in androgenic symptoms in females; acne, scalp hair loss, mood instability
- In males, excess 5 α -reductase activity and DHT can be associated with prostate concerns, scalp hair loss




Estrogens

Estrogen Neighborhood

- 🌱 Aromatase activity-
 - 🌱 E1 and E2
- 🌱 Evaluation of metabolites
 - 🌱 Phase 1 metabolism (-OH)
 - 🌱 Phase 2 metabolism (-M)
- 🌱 Evaluation of methylation potential (COMT enzyme)
- 🌱 Evaluation of potential for DNA mutagenesis
- 🌱 Evaluation of the metabolism of endogenous and exogenous hormones (BHRT)


2-OH: generally considered safest
 4-OH: potential for DNA damage
 16-OH: considered highly estrogenic



2-M-E1/2-OH-E1

Methylation Activity

4-M-E1/4-OH-E1




Aromatase / CYP19A1

CYP19A1

(-) Chrysin, Zn, nettles, damiana, grape seed, ECGC, resveratrol, licorice, flavonoids, anti-inflammatories

(+) Coleus forskohlii, inflammation, stress, Zn deficiency, excess adipose, high insulin, alcohol





Key Enzymes affecting Estrogen Metabolism


- 🔑 CYP1A1
- 🔑 CYP1B1
- 🔑 CYP3A4
- 🔑 COMT


Phase 1

🔗 **Phase 1 enzymes break down *fat-soluble molecules*.**

Phase 1 cytochrome P450 enzymes biotransform harmful substances by adding a hydroxyl, carboxyl, or amino group via either oxidation, reduction, or hydrolysis reactions (REDOX).

Phase 1 toxins - xenobiotics, antibiotics, ***steroid hormones***, pharmaceuticals.






ESTROGENS

The bar graph represents the relationship of the catechol estrogens (2-OH-E1, 4-OH-E1, 16-OH-E1) to each other. The expected percentage for each is represented by the shaded area.



The pathway illustrates phase 1 and phase 2 metabolism of both E1 and E2. Phase 1 metabolites, also known as catechol estrogens, are active and can induce estrogenic actions. Phase 2 metabolism gives insight into a patient's ability to methylate, or potentially inactivate harmful metabolites.



2-OH: generally considered safest

4-OH: potential for DNA damage

16-OH: considered highly estrogenic




Phase 1 Detox-Estrogen Metabolism

“Catechol Estrogens”

- 2-OH E1,E2- “safe” estrogen (CYP1A1)
- 4-OH E1,E2- DNA damage potential (CYP1B1)
- 16-OH E1- proliferative (CYP3A4)

CYP1A1 and CYP1B1



Mnemonic


🍌 CYP1A1 – A is awesome!!!

🍌 CYP1B1 – B is bad!!!



CYP1A1

- Healthy pathway of phase 1
 - Hydroxylates E1 to 2-OH E1
- Also, metabolizes polycyclic aromatic hydrocarbons (PAHs), such as benzo[a]pyrene, released from burning coal, oil, gasoline, trash, tobacco, wood, or charcoal-broiled meat.
- Mainly involved in detoxification, though can sometimes activate cancer-promoting substances such a PAHs PMID: 19531241 PMID: 28074113
- Can also activate aflatoxin B1 (cancer promoting) and tobacco-related carcinogens PMID: 19531241




3,3'-Diindolylmethane (DIM)

- DIM supplementation has been shown to enhance the 2-hydroxylation of estrogen, which is thought to contribute to anti-inflammatory effects in hormone-responsive cell lines
- DIM-treated subjects, relative to placebo, showed a significant **increase in levels of 2-OHE1**
- The influence of DIM results in **reduced production of the carcinogenic 4-hydroxyesterone (4OHE₁)**

- Szaefer H, Licznerska B, Krajka-Kuźniak V, Bartoszek A, Baer-Dubowska W. Modulation of CYP1A1, CYP1A2 and CYP1B1 expression by cabbage juices and indoles in human breast cell lines. *Nutr Cancer*. 2012 Aug;64(6):879-88. doi: 10.1080/01635581.2012.690928. Epub 2012 Jun 20. PMID: 22716309.
- Vivar OI, Saunier EF, Leitman DC, Firestone GL, Bjeldanes LF. Selective activation of estrogen receptor-beta target genes by 3,3'-diindolylmethane. *Endocrinology*. 2010 Apr;151(4):1662-7. doi: 10.1210/en.2009-1028. Epub 2010 Feb 16. PMID: 20160136; PMCID: PMC2850231.
- Dalessandri KM, Firestone GL, Fitch MD, Bradlow HL, Bjeldanes LF. Pilot study: effect of 3,3'-diindolylmethane supplements on urinary hormone metabolites in postmenopausal women with a history of early-stage breast cancer. *Nutr Cancer*. 2004;50(2):161-7

CYP1B1

- Takes estrogen down the 4 OH pathway
 - Pushes E1 to 4-OH E1 and E2 to 4-OH E2
- CYP1B1 is more readily found in tumor tissue compared to normal.
- Can activate cancer promoting compounds
- Accumulating evidence indicates that modulation of CYP1B1 can decrease adipogenesis and tumorigenesis, and prevent obesity, hypertension, atherosclerosis, and cancer. PMID: 28322972
- Given the role of CYP1B1 in pro-carcinogen and estrogen metabolism, polymorphisms in CYP1B1 could result in modifications in its enzyme activity and subsequently lead to hormone-mediated carcinogenesis.
- Increases fat uptake and can lead to factors of metabolic syndrome



4-OH Estrogens Increase Cancer Risk

- 4-Hydroxylated catechol estrogens possess carcinogenic potential due to their ability to cause DNA damage by forming depurinating adducts, which in turn, generate mutations with subsequent oxidative damage and initiation of breast cancer.
- In microsomal preparations of human mammary adenocarcinoma, formation of 4-**hydroxy**estradiol was four times higher than 2-**hydroxy**estradiol formation, indicating that the ratio of 4-/2-hydroxyestradiol may be used as a biomarker for detection of malignant breast tumors.
- On the other hand, it has been suggested that the 4-**methoxy**estrogens prevent oxidative metabolism of estradiol and oxidative DNA damage.


Cancer chemoprevention revisited: Cytochrome P450 family 1B1 as a target in the tumor and the microenvironment

- CYP1B1 converts E_2 and E_1 into 4-hydroxy estradiol (4-OHE₂) or 4-hydroxy estrone (4-OHE₁). 4-OHE₂/E₁ can be further converted into highly potent carcinogens, such as quinones and semiquinones, leading to reactive oxygen species and DNA adducts.
- In contrast, CYP1A1 converts E_2 and E_1 into 2-hydroxy estradiol (2-OHE₂) and 2-hydroxy estrone (2-OHE₁). These metabolites have low affinity for estrogen receptors; although they might also undergo metabolic redox cycling to generate free radicals and reactive semiquinone/quinone intermediates, they are rapidly converted by the enzyme COMT into 2-methoxy estradiol (2-ME₂) or 2-methoxy estrone (2-ME₁), which are known to be protective against tumor formation.
- CYP1B1 is therefore implicated in the development and progression of estrogen-mediated tumors, as its overexpression is responsible for estrogens conversion into toxic metabolites

CYP3A4

Moves E2 to E3 and E1 to 16OHE1
(and metabolizes testosterone)

- 🌱 CYP3A4 metabolizes about half of all drugs on the market
- 🌱 Mostly found in the liver, but most active in the gut
- 🌱 Metabolizes many internal compounds such as cholesterol, fatty acids, prostaglandins, leukotrienes, retinoids and biogenic amines PMID: 25332983
- 🌱 Detoxifies bile acids PMID: 25332983
- 🌱 Partially degrades vitamin D PMID: 22985909
- 🌱 The activity of this enzyme varies and is affected by health, smoking, diet, medication, hormones, and genetics



Estriol

- Estriol is the main estrogen in pregnancy but has received less attention outside gestation
- Non-proliferative estrogen
- Does not convert to E1 or E2
- Immunomodulatory benefits- autoimmune, inflammatory, and neurodegenerative conditions
- Treatment of menopausal symptoms, osteoporosis, cancer, hyperlipidemia, vascular disease, and multiple sclerosis.

16-OH E1-Proliferative / Tissue Growth


- Good for bones
- Not so good for breast/fibroids/endometriosis
- Possibly more estrogenic symptoms
- Binds most strongly to the estrogen receptor compared to 2OH and 4OH
 - Though much weaker than E1 and E2

CYP3A4

CYP3A4


(-) Grapefruit, aloe, polyphenols (i.e. resveratrol), flavonoids, ECGC, coffee, fennel, black pepper, licorice, chrysin, quercetin

(+) St. John's wort, capsaicin, valerian, ginkgo biloba, fatty acids, Vit D



How do you lower 16-OHE1?

- 🌿 DIM
- 🌿 Resveratrol
- 🌿 Quercetin
- 🌿 Sulforaphane
- 🌿 Rosemary




2, 4, 16-OH E1 Compared to Each Other

2-OH: generally considered safest

4-OH: potential for DNA damage

16-OH: considered highly estrogenic




Phase 2 Detoxification

🌱 **Phase 2 enzymes** break down compounds into *water-soluble compounds* for safer excretion.

- 🌱 Conjugated with a water-soluble or hydrophilic substance in phase 2.
- 🌱 This is completed by adding various groups such as sulfur, glutathione, glucuronic acid, amino acids, or methyl or acetyl groups.
- 🌱 Phase 2 enzymes are: UDP-glucuronosyl transferases, sulfotransferases, glutathione transferases, amino acid transferases, N-acetyl transferases, and methyltransferases/catechol-O-methyltransferase.


Phase 2

🔗 If the 2-OH and 4-OH estrogens are not methylated, they can be converted to highly reactive molecules that could potentially damage DNA.





COMT (Catechol-O-Methyltransferase)



COMT

(-) High sucrose diet, leptin resistance, inflammation, bisphenol and PCBs

(+) B12, folate, Mg, methionine, SAMe, betaine, TMG, resveratrol, citrus, rosemary, DIM, rooibos, curcumin

Estrogen Metabolism: Phase 2 Detox

2-M-E1/2-OH-E1

Methylation
Activity

4-M-E1/4-OH-E1

Oxidation or Methylation?

- 🌀 Methylation (COMT)
 - 🌀 Inactivates catechol estrogens
 - 🌀 Increases solubility/prepare intermediates for renal and biliary excretion
 - 🌀 methoxy metabolites → are stable and safely leave the body
- 🌀 Oxidation
 - 🌀 CYP peroxidase → quinone and semiquinone formation
 - 🌀 Glutathionization (GST)
 - 🌀 Glucuronidation (UGT)
 - 🌀 Sulfonation (SOD2)
 - 🌀 Acetylation

COMT

(-) High sucrose diet, leptin resistance, inflammation, bisphenol and PCBs


(+) B12, folate, Mg, methionine, SAmE, betaine, TMG, resveratrol, citrus, rosemary, DIM, rooibos, curcumin



2-M-E1:2-OH-E1 ⁺	(COMT/Methylation activity)	0.22		0.08 – 0.60
2-M-E2:2-OH-E2 ⁺	(COMT/Methylation activity)	0.13		0.06 – 0.80
4-M-E1:4-OH-E1 ⁺	(COMT/Methylation activity)	0.0087		0.004 – 0.10

Estrogen Quinones and Breast Cancer

- Un-methylated 2-OH and 4-OH metabolites can form quinones / semi-quinones, resulting in DNA adduct formation/DNA damage
- Estrogen quinones are highly carcinogenic
- Potentiated by oxidative damage
- Genotoxic: DNA breaks/adducts
- Mutagenic: Depurinated DNA
- Carcinogenic: high levels in BC and HR breast tissue
 - Induction and propagation of BC
 - Increased risk of hormone responsive BC
- “Increased amounts of estrogen-DNA adducts are found not only in people with several different types of cancer but also in women at high risk for breast cancer, indicating that the formation of adducts is on the pathway to cancer initiation.”**




Quinone reduction


NAC

- Reduce estrogen semiquinones back to catechol estrogens
- Primary effect is to react with quinones to form conjugates preventing the formation of estrogen-DNA adducts.

Resveratrol

- Reduce catechol estrogen semiquinones back to catechol estrogens
- Induce the estrogen-protective enzyme quinone reductase
- Modulates CYP1B1, thereby reducing its activity and thus the formation of 4-OHE1(E2)






Phase 3 Detox

Beta glucuronidase: Elimination or Re-circulation

ESTROGEN METABOLISM IN THE GUT



- Conjugated estrogens pass into the gut for elimination
- Beta-glucuronidase can deconjugate estrogens
- Unconjugated estrogens reabsorbed into circulation
- High levels of beta-glucuronidase could mean more recirculating of estrogens
- Beta-glucuronidase increases with: low fruits and veggies, low fiber, high sugar, processed foods, SAD diet, alcohol, toxicants, antibiotics



Order: 999999-9999
 Client #: 999999
 Doctor: Sample Doctor, ND
 Doctors Data Inc
 123 Main St.
 St. Charles, IL 60174 USA

Patient: Sample Patient
 Id: 999999 Client Ref: ICL220706054
 Age: 43 DOB: 10/25/1978
 Sex: Female

Sample Collection Date/Time
 Date Collected 07/04/2022
 Date Received 07/07/2022
 Date Reported 07/17/2022
 Specimens Collected 3

Microbiome Abundance and Diversity Summary

The abundance and diversity of gastrointestinal bacteria provide an indication of gastrointestinal health, and gut microbial imbalances can contribute to dysbiosis and other chronic disease states. The GI360™ Microbiome Profile is a gut microbiota DNA analysis tool that identifies and characterizes more than 45 targeted analytes across six Phyla using PCR and compares the patient results to a characterized normobiotic reference population. The web chart illustrates the degree to which an individual's microbiome profile deviates from normobiosis.




Dysbiosis and Diversity Index

These indexes are calculated from the results of the Microbiome Profile, with scores ranging from 1 to 5, and do not include consideration of dysbiotic and pathogenic bacteria, yeast, parasites and viruses that may be reported in subsequent sections of the GI360™ test.

The Dysbiosis Index (DI) is calculated strictly from the results of the Microbiome Profile, with scores from 1 to 5. A DI score above 2 indicates dysbiosis; a microbiota profile that differs from the defined normobiotic reference population. The higher the DI above 2, the more the sample deviates from the normobiotic profile. The dysbiosis test and DI does not include consideration of dysbiotic and pathogenic bacteria, yeast, parasites and viruses that may be reported in subsequent sections of the GI360™ test.

A diversity score of 3 indicates an expected amount of diversity, with 4 & 5 indicating an increased distribution of bacteria based on the number of different species and their abundance in the sample, calculated based on Shannon's diversity index. Scores of 1 or 2 indicate less diversity than the defined normobiotic reference population.



Key Findings

- Butyrate producing bacteria Clostridioides difficile (Toxin A/B), Detected
 - Gut barrier protective bacteria RBC, Abnormal
 - Gut intestinal health marker Candida albicans, Cultured
 - Pro-inflammatory bacteria
 - Gut barrier protective bacteria vs. opportunistic bacteria
- = Expected = Imbalanced

Stool testing: beta glucuronidase

Digestion / Absorption	Result	Unit	L	WRI	H	Reference Interval
Elastase	>500	µg/mL				> 200
Fat Stain	None					None – Few
Carbohydrates [†]	Negative					Negative
Inflammation	Result	Unit	L	WRI	H	Reference Interval
Lactoferrin	1.6	µg/mL				< 7.3
Lysozyme*	369	ng/mL				≤ 500
Calprotectin	<10	µg/g				≤ 50
Immunology	Result	Unit	L	WRI	H	Reference Interval
Secretory IgA*	34.8	mg/dL				30 – 275
Short Chain Fatty Acids	Result	Unit	L	WRI	H	Reference Interval
% Acetate [‡]	62					50 – 72
% Propionate [‡]	17					11 – 25
% Butyrate [‡]	17					11 – 32
% Valerate [‡]	3.7					0.8 – 5.0
Butyrate [‡]	2.4	mg/mL				0.8 – 4.0
Total SCFA's [‡]	14	mg/mL				5.0 – 16.0
Intestinal Health Markers	Result	Unit	L	WRI	H	Reference Interval
pH	6.1					5.8 – 7.0
β-glucuronidase*	1508	U/L				100 – 1200
Occult Blood	Negative					Negative

Consider treating phase 3 (elimination) first

- 🌿 Hydration PMID: 20646222
- 🌿 Fiber PMID: 26026145
- 🌿 Herbs to prevent/treat constipation (aloe, ginger, Avipattikar, Triphala) PMID: 30680163 PMID: 28696777
- 🌿 Magnesium
- 🌿 Movement
- 🌿 Calcium D-Glucarate PMID: 2346674

Assessing oxidative damage: 8-hydroxy-2-deoxyguanosine(8-OHdG)


- Not an estrogen metabolite, but pairs nicely with metabolite information
- Lipids of cellular membranes, proteins, and DNA can incur permanent oxidative damage.
- In nuclear and mitochondrial DNA, 8-OHdG is one of the predominant forms of free radical-induced oxidative lesions and widely used as a biomarker for oxidative stress and carcinogenesis.
- 8-OHdG may be helpful in identifying or confirming DNA damage from harmful metabolites (quinones and DNA adducts)

Oxidative Stress Metabolite	Result	Unit	L	WRI	H	Reference Interval
8-hydroxy-2'-deoxyguanosine [†]	9.08	ng/mg Creat/Day				0 – 7.5


Oxidative stress and disease states

- Chronic stress, Cortisol elevation
- Inflammation
- Insomnia in postmenopausal women
- Degenerative diseases (rheumatoid arthritis, Parkinson's disease, Huntington's disease, Alzheimer's disease)
- Chronic fatigue syndrome
- Major depression
- Hypertension, Cardiovascular disease
- Diabetes type II
- Cystic fibrosis
- Psoriasis
- Chronic hepatitis
- Gastritis
- Irritable bowel disease
- Pancreatitis
- Cancer
- Acute viral infection
- Copper implants
- Toxic exposures (tobacco smoke, methamphetamines, asbestos, heavy metals, polycyclic hydrocarbons).


THE TESTS




Important, clinically relevant markers are color-coded to respond to an individual's imbalances, providing an at-a-glance format on the summary page.




Patient #1



Patient #2



Patient #3






Hormone & Urinary Metabolites Assessment Profile

Order: 999999-9999
 Test: X999999-9999-1
 Client #: 999999
 Doctors Data Inc
 123 Main St.
 St. Charles, IL 60174 USA

Patient: Sample Patient
 Id: 999999
 Age: 61 DOB: 01/01/1960
 Sex: Female
 Body Mass Index (BMI): 25
 Menopausal Status: Post-menopausal

Sample Collection Date: 01/01/2020
 Dinnertime: 12:00
 Bedtime: 12:00
 Waking: 12:00
 2 Hr. Post Waking: 12:00
 Collection Period: Morning
 Date Received: 01/01/2020
 Date Reported: 01/01/2020



ESTROGENS

The bar graph represents the relationship of the catechol estrogens (2-OH-E1, 4-OH-E1, 16-OH-E1) to each other. The expected percentage for each is represented by the shaded area.

The pathway illustrates phase 1 and phase 2 metabolism of both E1 and E2. Phase 1 metabolites, also known as catechol estrogens, are active and can induce estrogenic actions. Phase 2 metabolism gives insight into a patient's ability to methylate, or potentially inactivate harmful metabolites.

2-OH: generally considered safest
 4-OH: potential for DNA damage
 16-OH: considered highly estrogenic

2-OH-E1 % EXPECTED	4-OH-E1 % EXPECTED
50-85	2-7

CORTICOIDS

11βHSD2 is responsible for the conversion of cortisol to cortisone. Inhibition of this enzyme may lead to the amount of cortisol being greater than cortisone, while increased enzyme activity can lead to higher levels of cortisone in comparison to cortisol.

KEY RELATIONSHIPS

The graphs to the right represent metabolism preference by key enzymes, indicated by the arrow.

Metabolites in the 5-alpha pathway are more androgenic than their 5-beta counterparts and can be responsible for androgenic symptoms even when hormone levels appear normal.

Aromatase is an enzyme found in the greatest amounts in peripheral fat tissue which can increase estrogens in both males and females.

4-OH-E1 is considered unfavorable due to its carcinogenic potential within breast and prostatic tissue as a reactive metabolite. When methylated by COMT, this reactive metabolite becomes stable and can be removed from the body.

5-α REDUCASE

5β preference vs. 5α preference

AROMATASE ACTIVITY

TESTOSTERONE vs. ESTRADIOL

ANDROSTENEDIONE vs. ESTRONE


COMT/METHYLATION ACTIVITY

4-HYDROXYESTRONE vs. 4-METHOXYESTRONE

ESTROGENS

The bar graph represents the relationship of the catechol estrogens (2-OH-E1, 4-OH-E1, 16-OH-E1) to each other. The expected percentage for each is represented by the shaded area.


The pathway illustrates phase 1 and phase 2 metabolism of both E1 and E2. Phase 1 metabolites, also known as catechol estrogens, are active and can induce estrogenic actions. Phase 2 metabolism gives insight into a patient's ability to methylate, or potentially inactivate harmful metabolites.



2-OH: generally considered safest

4-OH: potential for DNA damage

16-OH: considered highly estrogenic










Order: 999999-9999
Test: X999999-9999-1
Client #: 999999
Doctors Data Inc
123 Main St.
St. Charles, IL 60174 USA

Patient: Sample Patient
Id: 999999
Age: 61 DOB: 01/01/1960
Sex: Female
Body Mass Index (BMI): 25
Menopausal Status: Post-menopausal

Sample Collection Date/Time
Dinnertime 12/30/2022 19:20
Bedtime 12/30/2022 22:30
Waking 12/31/2022 07:00
2 Hr. Post Waking 12/31/2022 09:30
Collection Period Multipoint daily
Date Received 01/01/2023
Date Reported 01/02/2023

Cortisol and cortisone are plotted on the same graph to quickly identify their relationship to each other.



Analyzed by DOCTOR'S DATA, INC. • 3755 Illinois Avenue, St. Charles, IL 60174-2420 USA • LAB DIR: Eric Roth, MD • CLIA ID: 14D0646470




Order: 999999-9999
Test: X999999-9999-1
Client #: 999999
Doctors Data Inc
123 Main St.
St. Charles, IL 60174 USA

Patient: Sample Patient
Id: 999999
Age: 61 DOB: 01/01/1960
Sex: Female
Body Mass Index (BMI): 25
Menopausal Status: Post-menopausal

Sample Collection Date/Time
Dinnertime 12/30/2022 19:20
Bedtime 12/30/2022 22:30
Waking 12/31/2022 07:00
2 Hr. Post Waking 12/31/2022 09:30
Collection Period Multipoint daily
Date Received 01/01/2023
Date Reported 01/02/2023

Important relationships are highlighted here, representing metabolism preference by key enzymes.



KEY RELATIONSHIPS

The graphs to the right represent metabolism preference by key enzymes, indicated by the arrow.

Metabolites in the 5-alpha pathway are more androgenic than their 5-beta counterparts and can be responsible for androgenic symptoms even when hormone levels appear normal.


Aromatase is an enzyme found in the greatest amounts in peripheral fat tissue which can increase estrogens in both males and females.

4-OH-E1 is considered unfavorable due to its carcinogenic potential within breast and prostatic tissue as a reactive metabolite. When methylated by COMT, this reactive metabolite becomes stable and can be removed from the body.

5-A REDUCTASE ACTIVITY




AROMATASE ACTIVITY




COMT/METHYLATION ACTIVITY



Analyzed by DOCTOR'S DATA, INC. • 3755 Illinois Avenue, St. Charles, IL 60174-2420 USA • LAB DIR: Erio Roth, MD • CLIA ID: 14D0646470



- Dynamic color-coded overview of the entire hormone steroid cascade displaying actual results for your patient.
- Colors are used to easily identify highs and lows for each metabolite.



Can't see the forest for the trees...

- Utilizing urinary hormone testing requires the practitioner to take a step back and look at the pattern of hormone secretion
- Metabolism is complex and will be different for every individual
- Attempting to bring every analyte within range is not necessarily the goal of this type of testing

Hormone and Urinary Metabolites Assessment Profile

Clinical and Therapeutic Considerations



Hormone & Urinary Metabolites Assessment Profile

Order: SAMPLE REPORT
 Client #: 12345
 Doctor: Sample Doctor
 Doctor's Data, Inc.
 3755 Illinois Ave.
 St. Charles, IL 60174

Patient: Sample Patient
 Age: 35
 Sex: Female
 Menopausal Status: Pre-menopausal

Sample Collection Date/Time: 06/14/2023
 Dinnertime: 06/14/2023
 Bedtime: 06/14/2023
 Waking: 06/14/2023
 2 Hr. Post Waking: Multipoint
 Collection Period: 06/15/2023
 Date Received: 06/16/2023
 Date Reported: 06/16/2023



ESTROGENS
 The bar graph represents the relationship of the 2-OH generally considered best 4-OH potential for DNA damage to each other. The expected percentage for each is represented by the shaded area.

The pathway illustrates phase 1 and phase 2 metabolism of both E1 and E2. Phase 1 metabolites, also known as catechol estrogens, are active and can induce estrogen actions. Phase 2 metabolites may be inactive or potentially inactive hormonal metabolites.



CORTICOSTEROIDS
 17βHSD is responsible for the conversion of cortisol to cortisone. Inhibition of this enzyme may lead to the amount of cortisol being greater than cortisone, while increased enzyme activity can lead to higher levels of cortisone in comparison to cortisol.




KEY RELATIONSHIPS
 The graphs to the right represent metabolites preferred by any enzyme, indicated by the arrows.

Metabolites in the 5-alpha pathway are more androgenic than their 5-beta counterparts and can be observed for androgenic symptoms when hormone levels appear normal.

Aromatase is an enzyme found in the granulosa cells in peripheral fat tissue which can convert androgens to both males and females. Increases in androgens is both males and females.

17βHSD is considered a cofactor due to its synergistic potential within breast and prostate tissue as a reactive metabolite. When metabolized by COMT, this reactive metabolite is metabolized to an inactive metabolite and can be removed from the body.

doctorsdata.com



aka Aromatase

Decrease 4-OH estrogen

Phase II detox support

Support 2-OH estrogen

1 CYP19A1
 (-) Chrysin, Zn, nettles, damiana, grape seed, EGCG, resveratrol, licorice, flavonoids, anti-inflammatories
 (+) Coleus forskohlii, inflammation, stress, Zn deficiency, excess adipose, high insulin, alcohol

2 HSD17B
 (-) Licorice, quercetin, apigenin, phytoestrogen, flax, tamoxifen
 (+) Grape seed, propolis, DHEA, coleus forskohlii, rooibos, flavonoids, alcohol, inflammation, abdominal adiposity

3 CYP21A2
 (-) Vit D, resveratrol, curcumin, rooibos, apigenin, isoflavonoids, DHEA, omeprazole, valproic acid, PDE inhibitors
 (+) Coleus forskohlii

4 CYP1A1
 (-) Resveratrol, EGCG, berries, St. John's wort, lycopene, propolis, grapefruit
 (+) DIM, rosemary, fish oil, tea, coffee, hops

5 SULT2A1/STS
 (-) Quercetin, licorice, spironolactone, testosterone, clomiphene, inflammation
 (+) Genistein

6 5β Reductase
 (-) Licorice, budesonide
 (+) NAD

7 CYP1B1
 (-) Apiaceae family, grapefruit, resveratrol, rosemary, St John's wort, quercetin, flavonoids, curcumin, EGCG, rooibos
 (+) THC, UV exposure, PAHs, PCBs, diesel exhaust, inflammation, insulin resistance, leptin resistance

8 COMT
 (-) High sucrose diet, leptin resistance, inflammation, bisphenol and PCBs
 (+) B12, folate, Mg, methionine, SAmE, betaine, TMG, resveratrol, citrus, rosemary, DIM, rooibos, curcumin

9 5α Reductase
 (-) Saw palmetto, nettles, pygeum, EGCG, progesterone, Zn, berberine, polyunsaturated fatty acids
 (+) DHEA, insulin resistance, obesity, PCOS, essential HTN, high carb diet, sodium restriction

10 HSD3A
 (-) Coumestrol, taxifolin, mirtazapine
 (+) Sulforaphane, oral progesterone, coleus forskohlii

11 HSD11B1
 (-) Physical activity, rooibos, holy basil, curcumin, bitter melon, EGCG, progesterone, coffee, hyperthyroid, estrogen excess
 (+) Glucocorticoids, inflammation, hypothyroid, impaired glucose tolerance, insulin resistance, visceral adiposity

12 HSD3B1/2
 (-) Resveratrol, isoflavonoids, progestins, phytoestrogens, trilostane, troglitazone, ketoconazole, metformin, chronic EtOH use, CAH
 (+) Fenugreek, pregnenolone, PCOS, high insulin, hyperadrenalism, hyperthyroid, inflammation, coleus forskohlii

13 CYP11B1/2
 (-) Flavonoids, isoflavones, DHEA, azoles
 (+) Vit D, Korean red ginseng, coleus forskohlii, heme, dexamethasone, progestins, tizanidine

14 CYP3A4
 (-) Grapefruit, aloe, polyphenols (i.e. resveratrol), flavonoids, EGCG, coffee, fennel, black pepper, licorice, chrysin, quercetin
 (+) St. John's wort, capsaicin, valerian, ginkgo biloba, fatty acids, Vit D

15 HSD11B2
 (-) Intense exercise, grapefruit, progesterone, licorice, fluoxymesterone, azoles
 (+) Glucocorticoids, NAD

16 CYP17A1
 (-) Resveratrol, curcumin, licorice, apigenin, isoflavones, spironolactone, azoles, metformin, nicotine, dioxins
 (+) Vit D, DHEA, coleus forskohlii, hyperglycemia, stress, alcohol, antiepileptics (high dose), PCBs

		coleus forskohlii	curcumin	DIM	EGCG	grape seed	licorice	quercetin	resveratrol	rooibos	rosemary
1	CYP19A1	+			-	-	-		-		
2	HSD17B	+				+	-	-			
3	CYP21A2	+	-						-	-	
4	CYP1A1			+	-				-		+
5	SULT2A1/STS						-	-			
6	5β Reductase						-				
7	CYP1B1		-		-			-	-	-	-
8	COMT		+	+					+	+	+
9	5α Reductase				-						
10	HSD3A	+									
11	HSD11B1		-		-					-	
12	HSD3B1/2	+							-		
13	CYP11B1/2	+									
14	CYP3A4				-		-	-	-		
15	HSD11B2						-				
16	CYP17A1	+	-				-		-		



SCIENCE + INSIGHT

A clinical laboratory providing innovative, accurate specialty testing since 1972.

> Create an Account

➔ Clinician Sign In

Search

TEST MENU

HOW TO ORDER

CREATE ACCOUNT

NEWS & KNOWLEDGE

ABOUT DOCTOR'S DATA

CONTACT US

Test Menu

Test Menu / Endocrinology / **Hormone & Urinary Metabolites Assessment Profile**

CardioMetabolic

Endocrinology

Environmental Exposure and Detoxification

Gastrointestinal Health

Nutritional Status

Toxic and Essential Elements

Collection Type



Hormone & Urinary Metabolites Assessment Profile

The Hormone and Urinary Metabolites Assessment Profile (HuMap™) provides a comprehensive overview of steroid hormones, their metabolites, and the efficiency of the enzymes that metabolize these hormones. This non-invasive test requires only 4 or 5 separate urine collections. Because the breakdown of hormones relies so heavily on processes within the liver, this test can also elucidate areas of interest as it pertains to conjugation of each metabolite. Additionally, testing urinary hormone metabolites can contribute to further understanding of endogenous hormone secretion, supplemental hormone utilization, enzyme activity, oxidative stress, and insight into how your body metabolizes hormones. [Patient FAQs and Best Practices for HuMap](#) [\[➔ LEARN MORE\]](#)



Find Out More

- View Sample Report
- Collection Instructions
- Resource Guide
- Detailed Information
- Best Practices for Specimen Collection
- Compare Endocrinology Profiles
- HuMap Clinical and Therapeutic Considerations

Useful for:

- **Female**
- Menopausal symptoms
- Breast Health
- Endometriosis
- PCOS
- PMS/PMDD
- **Male**



QUESTIONS?

